

Sling Guide

Beaucare Slings

Our Sling Range

Our sling range focuses on comfort, durability and fit. This guide is designed to aid you to choose the most suitable sling for your needs and environment.

- Choose between 10 designs, 3 fabric choices and 6 sizes
- Colour-coded edges and loops
- Both loop and clip fixings available
- Embroidered labels that last the life of the sling



Our team can help you select the right sling

Vicki Crombleholme
Training and Compliance Manager



Scan here for **Sling Options** and **Prices**



Attachments & Fabrics

Attachment Types



Loops


Various coloured loop straps provide adjustment and correct positioning angle for the user.



Clips


Should your hoist attach using clip fixings we are able to have compatible slings made to order.

Fabrics



Poly

- Durable, low friction fabric
- Ideal for general transfers
- Withstands high laundry temperatures



Mesh

- For use when bathing & showering
- Quickly drains water, facilitating fast drying



Spacer

- Designed to be used and left in place between transfers
- Withstands high laundry temperatures



220 MAX. (35 STONES)

Do not exceed sling or hoist capacity


Sling Range

For bespoke requirements or bariatric sizing, please **contact your Account Manager or scan here**




Universal Sling

A standard quick-fit, universal sling that works for 85% of individuals. Offers a more general solution to individuals lifting needs.



Universal Deluxe Sling


Featuring a head support and padding on the leg area, this sling offers a higher level of comfort and support than a standard universal sling.



Classic Hammock Sling

Designed to provide support and comfort by reducing hip movement.


Optional - Head Support



Deluxe Hammock Sling

A hybrid sling designed to offer the full support of a hammock sling with the function of a universal. Apply and remove while seated.


Optional - Head Support



Hammock Sling

Hammock style long seat sling. Designed for amputees and individuals who are not comfortable in a divided leg sling.

Optional - Head Support



In Chair Hammock Sling

Provides a secure fit with safety pommel. Suitable for those with amputations.




In Chair Split Leg Hammock Sling

Split leg design allowing you to fold away the leg supports preventing build up of heat and adding to user comfort.




Toileting Sling

Designed with a larger aperture which allows for easy access during toileting. Clipped belt with a fastening at the waist.




Transporting Sling

Designed to offer greater support for an individual using a stand aid hoist. Clipped belt with a fastening at the waist.





Stand Aid

Designed to be used in conjunction with a stand aid hoist. Clipped belt with a fastening at the waist.




Size Guide



Use The Size Table Below For:

Toileting Sling
Transporting Sling
Stand Aid Sling

SLING SIZE	USER HEIGHT*	USER WIDTH*
X SMALL	145-155	70-85
SMALL	150-170	75-90
MEDIUM	160-180	85-100
LARGE	160-185	90-115
X LARGE	170-195	100-125




Use The Size Table Below For:

Universal Sling
Deluxe Universal Sling
Classic Hammock Sling
Deluxe Hammock Sling
In Chair Hammock Sling
In Chair Split Leg Hammock Sling

SLING SIZE	SLING WIDTH*	USER BACK*	USER WIDTH*
X SMALL	77	49-52	45-54
SMALL	86	56-58	54-60
MEDIUM	96	62-68	60-72
LARGE	116	67-73	72-82
X LARGE	129	70-73	82-100

*All measurements in cm




Use The Size Table Below For:

Hammock Sling

SLING SIZE	SLING WIDTH*	USER BACK*	USER WIDTH*
X SMALL	81	104	45-54
SMALL	86	109	54-60
MEDIUM	96	114	60-68
LARGE	106	122	68-76
X LARGE	116	129	76-92

*All measurements in cm



Hornbeam Square West, Hornbeam Park, Harrogate,
North Yorkshire, HG2 8PA

01423 873666 | www.beaucare.com | sales@beaucare.com